SILICONE EAR PLUGS

USER INSTRUCTIONS & CARE GUIDE

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This product complies to the General Product Safety Directive (2001/95/EC) and EN 352-2:2020.
This product is not considered as a personal protective equipment.

FITTING INSTRUCTIONS

- Fit one ear plug at a time.
- Before fitting, dry hands and ears.
- In your hands, roll the ear plug into a stretched oval shape.
- Hold the top of the ear and pull gently upwards and outwards to open the ear canal.
- Place the end of the ear plug into your ear.
- Press down and hold your finger or thumb over the ear plug for 10-15 seconds, until it flattens and moulds to form a snug seal.
- This fitting process will need to be repeated for each use.
- If the ear plug is too big, remove a small piece and repeat the fitting process.

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Distributed by:
Mendota Merchants LLC
8 The Green, Suite 7347
Dover, DE 19901

SILICONE EAR PLUGS

USER INSTRUCTIONS & CARE GUIDE

Storage advice: Store in cool, dry place.

This product complies to the General Product Safety Directive (2001/95/EC) and EN 352-2:2020.

This product is not considered as a personal protective equipment.

FITTING INSTRUCTIONS

- Fit one ear plug at a time.
- Before fitting, dry hands and ears.
- In your hands, roll the ear plug into a stretched oval shape.
- Hold the top of the ear and pull gently upwards and outwards to open the ear canal.
- Place the end of the ear plug into your ear.
- Press down and hold your finger or thumb over the ear plug for 10-15 seconds, until it flattens and moulds to form a snug seal.
- This fitting process will need to be repeated for each ear.
- If the ear plug is too big, remove a small piece and repeat the fitting process.

SAFETY WARNINGS & RECOMMENDATIONS

- If you experience excessive pressure on the ear drum, then remove the ear plug and repeat fitting process.
- Use ear plugs only as instructed. If the ear plug is made too small, the ear plug can become lodged in the ear and cause harm.
- Ear plugs should be kept out of the reach of small children. If being used by children, they should be supervised by an adult. Swallowing, or inserting into the nose, could interfere with breathing and result in serious injury.
- Avoid inserting the ear plugs too far into the ear because they could become difficult to remove.
- Dispose of the ear plugs if they start to break up or disintegrate.
- The ear plugs are reusable, so to keep them clean and increase their longevity, wash regularly with warm water and allow to dry before reusing.
- Although reusable, please replace ear plugs regularly to avoid irritation or infection.

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